



**inhealth™ Fiab ECG Foam Monitoring Electrode**  
**- Snap 48x50mm Pkt/50**



<b>ITEM DESCRIPTION:</b>	Inhealth Fiab ECG Foam Monitoring Electrode - Snap
<b>SIZE:</b>	48 x 50mm
<b>PACK UOM:</b>	Box/50pcs
<b>PRODUCT CODE/S:</b>	225309
<b>MATERIAL:</b>	Adhesive Foam Backing

**inhealth™ Fiab ECG Non-Woven Electrode**  
**- Snap 55mm Pkt/50**



<b>ITEM DESCRIPTION:</b>	Inhealth Fiab ECG Non-Woven Electrode – Snap
<b>SIZE:</b>	55mm
<b>PACK UOM:</b>	Box/50pcs
<b>PRODUCT CODE/S:</b>	225310
<b>MATERIAL:</b>	Adhesive non-woven backing

FIAB electrodes provide the highest quality and performance during ECG tests. stress test, long and short monitoring, and Holter test.

- Solid gel has a special formula that improves adhesion and easy application. Adhesive backing is made in foam or non-woven.
- Conductive adhesive provides an excellent ECG reading and the soft medical foam or the non-woven with a sensitive adhesive ensure a secure comfortable contact.
- Connection is a standard snap adapter.
- For use on adults.

